

The Foot Strap from Morrison Medical

The Foot Strap was designed to prevent patients from sliding on a backboard while being transported in an inclined or declined position, commonly seen in stairways and confined areas. This product will allow your patient to feel and be secure while maintaining the proper spine and neck immobilization that has been applied.



- Black strip on inside of strap acts as guide for quick and easy application.
- Place black strip across arches of feet; criss-cross strap over the ankles; loop ends through hand holes of backboard and fasten buckles; pull to tighten.

- 100% x-ray translucent.
- Single use, disposable.
- 2" polypropylene webbing.

Item #1000: The Foot Strap

1 piece disposable poly with plastic buckle, each

*For more information, call 800-438-6677,
or visit www.MorrisonMed.com.*