


Patho-Shield® Straps

Morrison Medical's Patho-Shield® Straps will outlive and outperform your current straps with less maintenance required. With no exposed fibers, these straps are 100% impervious to blood and other bodily fluids, making them quicker and easier to clean than standard plastic and cloth straps. The longevity of these straps means fewer replacements over time, saving you money. Patho-Shield® Straps are extremely strong – with a break strength of more than 1,800 pounds – yet adjust easily and stay very flexible in extreme temperatures. Make your Patho-Shield® straps identifiable by having us print your company name or logo directly on the strap. Don't see the strap combination you're looking for? Call us with your specifications and we'll do our best to accommodate your needs.



NEW COLOR NOW AVAILABLE
SAFETY GREEN



Patho-Shield® Straps are available in orange, yellow, black, and NEW SAFETY GREEN.

Please Call for Pricing

- #1300P: 2 piece, metal push button buckle with loop roller buckle ends, 5'
 - #1301P: 2 piece, metal push button buckle with loop roller buckle ends, 7'
 - #1308P: 1 piece, metal push button buckle, 7'
 - #1310P: 2 piece, metal push button buckle with swivel speed clip ends, 5'
 - #1311P: 2 piece, metal push button buckle with swivel speed clip ends, 7'
 - #1360P: 2 piece, plastic side release buckle with swivel speed clip ends, 5'
 - #1385P: Shoulder harness system with loop roller buckle ends
 - #1317P: 3' extension strap
 - #1318P: 4' extension strap
- Additional configurations available.*



Shoulder harnesses available.



How do Morrison Medical's Patho-Shield® Straps compare to the Biosafe® restraints ?

	Morrison Medical's Patho-Shield® Straps	vs.	Biosafe® restraints
Break Strength	Over 1,800 lbs.	Superior ★	1,200 lbs.
Flexibility	Remain flexible in extreme high and low temperatures.	★	Become stiff and difficult to work with in low temperatures.
Adjustability	Quickly and easily adjusts.	★	Difficult to pull to adjust, creating discomfort for patient and delaying transport.

For more information, call 800-438-6677, or visit www.MorrisonMed.com.